


Distortion-free, multi-shot EPI with Blip-Up and Down Acquisition (BUDA) and structured low-rank reconstruction

Contact: Congyu Liao, [cyliao@stanford.edu](mailto:cyliao@stanford.edu)

# Blip-up and down acquisition/reconstruction framework

Blip-up Recon Blip-down Recon Distortion-free Recon




Joint Reconstruction

$$\min_x \sum_t \|F_t E_t C x_t - d_t\|_2^2 + \lambda \|\mathcal{H}(x)\|_*$$

parallel imaging with field map  $E$

low-rank constraint

Iterative reconstruction




d1: shot1

blip-up

d2: shot2

blip-down




low-rank on Hankel matrix

(A)

(B)

# Acquired in-vivo datasets



Download these datasets from :

<https://drive.google.com/file/d/1DAeTW2iLgQx8tjZT3QqZ7F3vvt6Vwt2Q/view?usp=sharing>

And then copy these datasets to `./data/` folder

# The scan protocol of BUDA data


- FOV: 220 mm \*220mm, matrix size: 176\*176
- In-plane resolution: 1.25 mm \*1.25mm
- Slice thickness: 2mm
- Number of slices: 57
- Number of coils: 32
- TR/TE = 2800/77 ms
- Effective echo-spacing = 0.37 ms
- Multi-band factor: 3 with CAPI-shift FOV/2 per shot
- In-plane acceleration factor :2 per shot
- partial Fourier factor: 0.75
- BUDA acquisition with AP shot kyshift 1 and PA shot kyshift 0
- Two shells: b= 1000 with 32 diffusion directions and b=2500 with 64 diffusion directions

# MATLAB offline reconstruction

- The reconstruction scripts have been tested on MATLAB 2015b.
- The reconstruction scripts have been tested on Linux OS.
- The scripts require BART toolbox and FSL software
  - BART toolbox 0.5.00 is recommended (<https://zenodo.org/record/3376744>)
  - FSL 6.0.4 is recommended (<https://fsl.fmrib.ox.ac.uk/fsl/fslwiki/FslInstallation>)
- **Warning: your workstation should have more than 64GB RAM and 50 GB storage to reconstruct the raw data!**

# Reconstruction scripts

- Step1: run **S1\_SMS\_SENSE.m** to reconstruct each shot separately
- Step2: run **S2\_TOPUP.m** to estimate the field maps between AP/PA shots
- Step3: run **S3\_SMS\_BUDA.m** to jointly reconstruct distortion-free BUDA-EPI
- Step4: run **S4\_LoadALLVolumes.m** to load all volumes and combines AP/PA shots after removing background phase of each shot



# S1\_SMS\_SENSE.m

- 1. file\_name

```
28 %data path
29 file_path = './data/';
30 file_name = 'meas_MID00913_FID12373_BUDA_1p25inplane_R2MB3s2_ky10_blk_34dir.dat'; % raw data of BUDA-EPI
31 evp_name = 'meas_MID01329_FID12797_evp_baseline_for_PAT2MB3_with_ICE_on.dat'; % scan protocols
32 gre_name = 'meas_MID00909_FID12369_gre_sens_2mm57slc.dat'; % gre scan for sensitivity map
33
```

- 2. The “CoilSense\_ESPIRIT3d” function is using BART toolbox. Make sure you installed BART successfully.

```
72 % load GRE data and calculate sensitivity maps
73 if (~exist(strcat(save_path,'sens_gre.mat'),'file'))
74 fGRE=[file_path, gre_name];
75 [Img_gre, whmcc]= load_GRE_data_prewht(fGRE,N,Ncc); % load ref data, do coil compression and noise prewhitening
76 % estimate ESPIRiT coil sensitivity
77 sens_gre = sq_CoilSense_ESPIRIT3d(permute(Img_gre, [1,2,4,3])); % calculate coil sensitivity using ESPIRiT
78
```

- 3. Open this function, make sure the environment setting and folder path are correct


- 4. The script will create SMS-SENSE reconstructed AP/PA shots and ghost-corrected undersampled k-space data

k\_ap\_pa\_dif8rf1.mat  
k\_ap\_pa\_dif9rf1.mat  
prot.mat  
sens\_gre.mat  
img\_ap\_pa\_dif10rf1.nii  
img\_ap\_pa\_dif11rf1.nii  
img\_ap\_pa\_dif12rf1.nii

```
1 function sens = CoilSense_ESPIRIT3d(img_3d)
2 % this function is for 3d data, and estimates sensitivities in 3d
3
4 % img_3d: x,y,slice,chan
5 % data_path : path you want to write data
6 % sens : final coil sens output including all slices
7
8
9 ImSize = size(img_3d);
10 % code only work if matrix is even in size so mod matrix if it is odd
11
12 if rem(ImSize(1),2) == 1
13     img_3d = cat(1,img_3d,zeros(1,ImSize(2),ImSize(3),ImSize(4)));
14 end
15
16 if rem(ImSize(2),2) == 1
17     ImSizeCurrent = size(img_3d);
18     img_3d = cat(2,img_3d,zeros(ImSizeCurrent(1),1,ImSizeCurrent(3),ImSizeCurrent(4)));
19 end
20
21 if rem(ImSize(3),2) == 1
22     ImSizeCurrent = size(img_3d);
23     img_3d = cat(3,img_3d,zeros(ImSizeCurrent(1),ImSizeCurrent(2),1,ImSizeCurrent(4)));
24 end
25
26
27 addpath('/usr/local/app/bart/bart-0.5.00/matlab')
28
29 setenv('TOOLBOX_PATH','/usr/local/app/bart/bart-0.5.00')
30
31 setenv('TEMP_PATH','/dev/shm/');
32
33
```

# S2\_TOPUP.m

- Make sure the FSL functions can be called on MATLAB
- This script will generate field maps as well as unwrapped EPI shots





# S3\_SMS\_BUDA.m

- BUDA settings:



```
38 % BUDA settings
39 - opt.step_size = 0.8;
40 - opt.num_iter = 150; % 100 iterations
41 - opt.tol = 0.25;
42 - opt.winSize = [7,7]; % kernel size
43 - opt.lambda_msl = 1.00; % lambda of henkel low rank matrix
44 - opt.esp=esp; % echo-spacing
45 - opt.PhaseShiftBase = pi; % CAIPI shift of CAIPI
46 - opt.show_mercy = show_mercy ; % show how many cores you wanna use
47 - opt.fista=1 ; % --1 use FISTA algorithm --0 use POCS algorithm
48 - opt.t_k=1;
49
```

- You can also use non-parallel computing version to debug BUDA code

```
69 - tic
70 - % BUDA recon w/o parallel computing --it's slow, you can see the iterations
71 - % [img_msl]=recon_SMS_BUDA_XCLv2_fista (k_ap, k_pa,ky_idx_ap,ky_idx_pa,sens_gre, AccY, AccZ,img_fieldmap,opt);
72 -
73 - % BUDA_recon with parallel computing
74 - [img_msl]=recon_SMS_BUDA_XCLv2_parfor_fista (k_ap, k_pa,ky_idx_ap,ky_idx_pa,sens_gre, AccY, AccZ,img_fieldmap,opt);
75 -
76 - toc
```


- This script will generate BUDA reconstructed images



```
img_buda_rf1_dif11.mat
img_buda_rf1_dif12.mat
img_buda_rf1_dif13.mat
img_buda_rf1_dif14.mat
img_buda_rf1_dif15.mat
img_buda_rf1_dif16.mat
img_buda_rf1_dif17.mat
img_buda_rf1_dif18.mat
img_buda_rf1_dif19.mat
img_buda_rf1_dif2.mat
img_buda_rf1_dif20.mat
img_buda_rf1_dif21.mat
```

# S4\_LoadALLVolumes.m

- load whole-brain volumes and do real-valued diffusion processing to combine BUDA shots
- This script will generate final BUDA results in ./recon/ folder.



# Reference

If possible, please include the following citations if you find our reconstruction useful:

- [1] Liao C, Bilgic B, Tian Q, Stockmann JP, Cao X, Fan Q, Iyer SS, Wang F, Ngamsombat C, Lo WC, Manhard MK. Distortion-free, high-isotropic-resolution diffusion MRI with gSlider BUDA-EPI and multicoil dynamic B0 shimming. Magnetic Resonance in Medicine. 2021
- [2] Liao, C., Cao, X., Cho, J., Zhang, Z., Setsompop, K. and Bilgic, B., Highly efficient MRI through multi-shot echo planar imaging. In Wavelets and Sparsity XVIII (Vol. 11138, p. 1113818). International Society for Optics and Photonics. 2019